

## PHARMACEUTICAL ANALYSIS REPORT

5 July 2021

DISCLAIMER: The results herein are gathered for primary research purposes at the University of Notre Dame. This report and any attachments are intended for internal evaluation between the University of Notre Dame and the Liberia Medicines & Health Products Regulatory Agency. The results are not intended to discredit or replace the certificate of analysis provided by the manufacturer of the finished pharmaceutical product. If the results need to be used in an official capacity, further testing should be conducted at a laboratory with the appropriate jurisdiction.

Specific questions should be directed to:

Prof. Marya Lieberman, PhD  
Department of Chemistry and Biochemistry  
University of Notre Dame  
Notre Dame, IN 46556  
USA  
mlieberm@nd.edu  
(574) 631-4665

Report Prepared by:

Name	Institution	Position
Marya Lieberman	University of Notre Dame	Principal Investigator
Kathleen Hayes	University of Notre Dame	DPAL Coordinator Analyst

## I. SUMMARY

According to the United States Pharmacopeia (USP), finished pharmaceutical dosage forms of doxycycline must contain 90.0%-120.0% of the stated amount of active pharmaceutical ingredient (API). The University of Notre Dame has analyzed 13 samples of doxycycline representing 5 manufacturers and 5 batch numbers. Each doxycycline dosage form analyzed met the requirement for API content.

<b>Active Ingredient</b>	<b>Number of Samples Analyzed</b>	<b>Number of Failing Samples</b>	<b>Comment on Medicine Quality</b>	<b>Average API Content (%)</b>
Doxycycline	13	0	All passed	97.3%

## II. RESULTS

### A. HPLC SCREENING

Thirteen doxycycline samples collected in 2020 by covert shoppers in Liberia were selected for HPLC assay. One tablet from each package was assayed. If the sample failed assay, it was retested. All thirteen samples of doxycycline passed assay. We analyzed doxycycline made by five manufacturers, representing 5 batch numbers. All the samples had been stored in our 4°C cold room.

Sample	P1 API Content (%)	P1 Retest API Content (%)	Brand	Manufacturer	Lot Number	Expiration Date
21L-0005	95.3	NA	DoxyNAN	Saintsun Pharma Co. Ltd, China	180944	9/1/2021
21L-0009	103.4	NA	Philco Doxy 100	Philco Pharma, Germany	1705044	9/1/2021
21L-0004	102.2	NA	DoxyLab	Laborate Pharmaceutical, India, Ltd	DAOC-003	10/1/2022
21L-0007	95.9	NA	DoxyLab	Laborate Pharmaceutical, India, Ltd	DAOC-003	10/1/2022
21L-0010	101.6	NA	DoxyLab	Laborate Pharmaceutical, India, Ltd	DAOC-003	10/1/2022
21L-0012	96.1	NA	DoxyLab	Laborate Pharmaceutical, India, Ltd	DAOC-003	10/1/2022
21L-0001	99.7	NA	DoxyNext	Osmed Formulation Pvt,Ltd, India	K40901	10/1/2022
21L-0003	121.6	105.3	DoxyNext	Osmed Formulation Pvt,Ltd, India	K40901	10/1/2022
21L-0008	91.5	NA	DoxyNext	Osmed Formulation Pvt,Ltd, India	K40901	10/1/2022
21L-0011	92.5	NA	DoxyNext	Osmed Formulation Pvt,Ltd, India	K40901	10/1/2022
21L-0013	96.2	NA	DoxyNext	Osmed Formulation Pvt,Ltd, India	K40901	10/1/2022
21L-0002	95.4	NA	Doxycycline	Soham Healthcare Inc, India	SC-027	4/1/2022
21L-0006	90.3	NA	Doxycycline	Sohan Healthcare Inc, India	SC-027	4/1/2022

### III. METHODOLOGY

The procedures for quantifying the active ingredients in doxycycline finished pharmaceutical pills was based on the monograph published in the United States Pharmacopeia (USP). At times, minor modifications were made to the assay methodology. However, the accuracy, precision, and linearity of the methodology was verified according to USP <1226>. Sampling for the Distributed Pharmaceutical Analysis Lab is performed on a per pill basis, not as an average of 10 pills, since samples are generally composed of one dosage unit of the pharmaceutical product.

#### A. HPLC SYSTEM SUITABILITY

External calibration standards were created from analytical grade reagents. Five injections of the external standard were required to show a peak area within 2% relative standard deviation (RSD), and the range of retention times had to be within 0.3 minutes. After every five samples, the standard was injected as a quality check. It must have assayed within 2% RSD of the 5 initial injections and been within the 0.3 minute time range. If a quality check failed, data after the last passed quality check was discarded. Additional system suitability requirements information is available upon request.

#### B. HPLC INSTRUMENT CONDITIONS

Instrument: Waters 2695 High Performance Liquid Chromatograph

Column: XTerra C8, 5  $\mu$ m, 4.6 mm x150 mm

Temperature: RT

Detector: Waters 2487, Analytical Wavelength: 360 nm

#### C. HPLC SAMPLE PREPARATION

Sample Concentration: 0.5 mg/mL

Sample Solvent: 18 MOhm water

Sample Injection Volume: 20  $\mu$ L

#### D. HPLC STANDARD CONCENTRATION

Standard Concentration: 0.5 mg/mL

Standard Solvent: 18 MOhm water

Standard: Doxycycline hyclate, Sigma Aldrich, Lot # 069M4014V

#### E. HPLC MOBILE PHASE

Isocratic

30% Acetonitrile

70% 0.1% Trifluoroacetic acid in 18M Ohm water (pH = 1.8)

## APPENDIX

The appendix includes all the raw HPLC data and spreadsheet calculations for the samples and associated calibration checks.

## APPENDIX

Data Collection Date: 30 June 2021

Data Collection Location: University of Notre Dame

Workbook	Date
9/30/AP067	9/30/2021

STD #	Doxycycline Peak Area	STD Conc (mg/mL)	Manufacturer	Lot #	Expiry date**	Mass Dose (g)	Volume (mL)	Used (mL)	Final Vol (mL)	Purity*	Conc (mg/mL)
1	13546	Doxycycline	Sigma Aldrich	069M4014V	current lot	0.01266	6.0	0.75	3.0	0.956	0.5042
2	13637										
3	13644										
4	13646										
5	13638										
AVG	13633	*Express as pure doxycycline									
SD	54	**If using USP std, an expiry date is not assigned until they move to the next batch of standard, at which time the previous batch gets a one year expiry date. Enter "current lot."									
%RSD	0.39										
Pass?	pass										

Sample	Dox. Peak Area	Pill Total Mass (g)	Sample Mass (g)	Volume (mL)	Dox Dosage (mg)	mg Dox	%Error Dox	Dox Pass?	Comment	Total Mass (g)	MT Cap Mass (g)
21L-0001	26 13728	0.30343	0.01260	8.0	100	100	-0.3	pass			
21L-0002	27 12582	0.25899	0.01235	8.0	100	95	-4.6	pass			
21L-0003	28 16749	0.32218	0.01318	8.0	100	122	21.6	fail			
21L-0004	29 17643	0.29337	0.01499	8.0	100	102	2.2	pass			
21L-0005	30 17655	0.23841	0.01263	8.0	100	95	-4.7	pass			

Cal. Check_1	Dox. Peak Area
25	13475
Error (%)	1.16
Pass?	pass

Sample	Dox. Peak Area	Pill Total Mass (g)	Sample Mass (g)	Volume (mL)	Dox Dosage (mg)	mg Dox	%Error Dox	Dox Pass?	Comment	Total Mass (g)	MT Cap Mass (g)
21L-0006	31 13635	0.27426	0.01076	8.0	100	90	-9.7	pass			
21L-0007	32 13995	0.29289	0.01265	8.0	100	96	-4.1	pass			
21L-0008	33 14461	0.28915	0.01352	8.0	100	92	-8.5	pass	0.35143	0.06228	
21L-0009	34 15178	0.19571	0.01275	12.0	100	103	3.4	pass			
21L-0010	35 14603	0.29394	0.01246	8.0	100	102	1.6	pass			

Cal. Check_2	Dox. Peak Area
1	13475
Error (%)	1.14
Pass?	pass

Sample	Dox. Peak Area	Pill Total Mass (g)	Sample Mass (g)	Volume (mL)	Dox Dosage (mg)	mg Dox	%Error Dox	Dox Pass?	Comment	Total Mass (g)	MT Cap Mass (g)
21L-0011	36 13449.000	0.31104	0.01338	8.0	100	93	-7.5	pass			
21L-0012	37 15599.000	0.29471	0.01446	8.0	100	96	-3.9	pass			
21L-0013	38 15480.000	0.30357	0.01457	8.0	100	96	-3.8	pass			
21L-0003	39 15107.000	0.32218	0.01368	8.0	100	105	5.3	pass	retest: pass		
NA	NA	NA	NA	NA	NA	NA	NA	NA			


  

Cal. Check_3	Dox. Peak Area
1	13418
Error (%)	1.64
Pass?	pass

\*USP requirements for tablets and capsules containing doxycycline: 90.0-120.0% of labeled amount on package  
 All yellow cells must have an entry, such as NA, to make the LC data explicitly clear.  
 This sheet only works for single-point calibration  
 Password to unlock sheet: password




Lab name: Notre Dame  
Method: Syringe Injection  
Data file: doxy std 01.CHR ()



Component	Retention	Area	Height
-	2.170	51.8948	4.688
-	2.485	728.0415	87.571
doxycycline	3.663	13548.6027	1140.246
-	4.170	247.9528	21.163
-	4.540	12.9701	1.385

14589.3517


Lab name: Notre Dame  
Method: Syringe Injection  
Data file: doxy std 02.CHR ()



Component	Retention	Area	Height
-	2.143	104.2030	10.326
-	2.460	747.8558	91.309
doxycycline	3.510	13697.7097	1189.476
-	4.085	243.0242	21.706
-	4.433	12.3834	1.390
14804.9859			




Lab name: Notre Dame  
Method: Syringe Injection  
Data file: doxy std 03.CHR ()



Component	Retention	Area	Height
-	0.043	10.7487	0.388
-	2.126	107.7888	10.875
-	2.443	723.2133	89.373
doxycycline	3.563	13644.8040	1207.149
-	4.045	235.8350	21.538
-	4.380	11.8878	1.383
14734.2776			


Lab name: Notre Dame  
Method: Syringe Injection  
Data file: doxy std 04 CHR ()



Component	Retention	Area	Height
-	2.143	95.8854	9.430
-	2.463	730.2420	91.338
doxycycline	3.590	13640.7316	1220.052
-	4.070	238.6380	21.656
-	4.403	12.4024	1.399
14717.8794			

Lab name: Notre Dame  
 Method: Syringe Injection  
 Data file: doxy std 05.CHR ( )


14717.1800



Component	Retention	Area	Height
-	2.136	85.6152	8.248
-	2.460	740.1230	93.381
doxycycline	3.586	13638.6108	1232.039
-	4.063	240.5496	21.954
-	4.403	12.3814	1.445


14717.1800

Lab name: Notre Dame  
 Method: Syringe Injection  
 Data file: 21L-0001.CHR ()




Component	Retention	Area	Height
-	2.143	58.6604	6.331
-	2.470	628.8977	80.956
doxycycline	3.593	13728.7828	1267.600
-	4.066	227.8032	20.338
1464.4, 1641			

Lab name: Notre Dame  
Method: Syringe Injection  
Data file: 21L-0002.CHR ()



Component	Retention	Area	Height
-	2.463	457.9970	60.442
-	3.120	566.4417	20.573
doxycycline	3.543	15592.0054	1519.924
-	4.033	406.7240	26.021
17023.1681			


Lab name: Notre Dame  
Method: Syringe Injection  
Data file: 21L-0003.CHR ()



Component	Retention	Area	Height
-	2.136	46.6402	5.034
-	2.460	541.6352	73.073
doxycycline	3.546	16749.9076	1580.234
-	4.040	273.1652	24.939


17811.2542

Lab name: Notre Dame  
 Method: Syringe Injection  
 Data file: 21L-0004.CHR ()



Component	Retention	Area	Height
-	2.120	55.6614	6.304
-	2.440	753.0838	99.843
doxycycline	3.536	17643.6508	1611.687
-	4.026	431.9266	28.670
1888.2523			


Lab name: Notre Dame  
Method: Syringe Injection  
Data file: 21L-0005.CHR ()



Component	Retention	Area	Height
-	2.210	54.1808	5.753
-	2.533	674.8254	90.528
doxycycline	3.516	17660.6199	1654.709
-	4.110	276.2508	26.019
-	4.433	15.1221	1.773
18670.9983			




Lab name: Notre Dame  
Method: Syringe Injection  
Data file: cal check 01.CHR ()



Component	Retention	Area	Height
-	0.048	12.3672	0.273
-	2.113	100.8480	8.956
-	2.443	754.8596	98.046
doxycycline	3.533	13475.3748	1300.686
-	4.013	232.9754	22.616
-	4.333	13.4082	1.628
14589.8191			

Lab name: Notre Dame  
Method: Syringe Injection  
Data file: 21L-0006.CHR ()


Report generated by: C:\Program Files\Agilent\ChemStation\Bin\GENSTAT.EXE



Component	Retention	Area	Height
-	2.110	44.7920	4.504
-	2.426	974.5054	129.598
-	3.090	1030.7985	36.577
doxycycline	3.503	18658.8510	1779.641
-	4.006	473.8195	30.757


2118.7644

Lab name: Notre Dame  
Method: Syringe Injection  
Data file: 21L-0007.CHR ()




Component	Retention	Area	Height
-	0.063	22.0513	0.684
-	2.123	79.0724	9.130
-	2.450	738.1389	99.161
doxycycline	3.546	13995.5390	1332.890
-	4.023	285.4408	20.729
15120.2424			

Lab name: Notre Dame  
Method: Syringe Injection  
Data file: 21L-0008.CHR ()




Lab name: Notre Dame  
 Method: Syringe Injection  
 Data file: 21L-0009.CHR ()



Component	Retention	Area	Height
-	2.100	94.4298	11.008
-	2.423	789.8538	108.344
doxycycline	3.486	15178.6525	1496.693
-	3.966	194.6700	18.590
	1625.5161		


Lab name: Notre Dame  
Method: Syringe Injection  
Data file: 21L-0010.CHR ()



202007.1020050000




Lab name: Notre Dame  
Method: Syringe Injection  
Data file: 21L-0011.CHR ()







Lab name: Notre Dame  
Method: Syringe Injection  
Data file: 21L-0012.CHR ()




Lab name: Notre Dame  
Method: Syringe Injection  
Data file: 21L-0013.CHR ()



Lab name: Notre Dame  
Method: Syringe Injection  
Data file: 21L-0003 try 2.CHR ()



Lab name: Notre Dame  
Method: Syringe Injection  
Data file: cal check 03.CHR ()



Component	Retention	Area	Height
-	2.040	115.2795	11.885
-	2.426	774.4180	104.105
doxycycline	3.500	13410.4159	1330.954
-	3.966	228.6334	22.538
-	4.270	12.2351	1.532
-	4.670	19.8840	1.519
14560.8839			